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
		A25	A15
GENERAL	Test name	AMYLASE DIRECT	AMYLASE DIRECT
	Analysis mode	kinetic mon.	kinetic mon.
	Sample type	serum	serum
	Units	U/L	U/L
	Reaction type	increasing	increasing
	Decimals	0	0
	Replicates	1	1
Name of assoc. constituent	-	-	
PROCEDURE	Type of reading	monoch.	monoch.
	Sample	6	6
	Reagent 1	300	300
	Reagent 2	-	-
	Washing	1.2	1.2
	Predilution factor	-	-
	Main	405	405
	Reference	-	-
	Reading 1	60 s	72 s
	Reading 2	195 s	216 s
Reagent 2	-	-	
Postdilution factor	2	2	
CALIBRATION	Type of calibration	multiple	multiple
	Calibrator replicates	3	3
	Blank replicates	3	3
	Calibration curve	-	-
OPTIONS	Blank absorbance limit	0.200	0.200
	Kinetic blank limit	-	-
	Linearity limit	1300	1300

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- Winn-Deen ES, David H, Sigler G and Chavez R. Development of a direct assay for α -amylase. *Clin Chem* 1988; 34: 2005-2008.
- Gella FJ, Gubern G, Vidal R, Canalias F. Determination of total and pancreatic α -amylase in human serum with 2-chloro-4-nitrophenyl- α -D-maltotriósido as substrate. *Clin Chim Acta* 1997; 259: 147-160.
- Gubern G, Balsells D, Ferragut R, Galán A, Gella FJ, et al. Procedimiento recomendado para la determinación en rutina de la concentración catalítica de α -amilasa en suero sanguíneo humano. *Quim Clin* 1996; 15: 51-52.
- Balsells D, Gella FJ, Gubern G, Canalias F. Reference values for α -amylase in human serum and urine using 2-chloro-4-nitrophenyl- α -D-maltotriósido as substrate. *Clin Chim Acta* 1998; 274: 213-217.
- Tietz Textbook of Clinical Chemistry, 3rd edition. Burtis CA, Ashwood ER. WB Saunders Co., 1999.
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Press, 1995.
- Tietz Textbook of Clinical Chemistry, 2nd edition. Burtis CA, Ashwood ER. WB Saunders Co., 1994.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.