# BIOCHEMISTRY CONTROL SERUM I, II

## SUERO CONTROL DE BIOQUÍMICA I, II

NIVEL I		NIVEL II	
COD 18005 20 x 5 mL	COD 18009 5 x 5 mL	COD 18007 20 x 5 mL	COD 18010 5 x 5 mL

### USO PREVISTO

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### INTENDED USE

The Biochemistry Control Serum is intended for quality control purposes only and is supplied with suggested intervals of acceptable values. For *in vitro* professional use only in the clinical laboratory.

### CONTENTS AND COMPOSITION

**Biochemistry Control Serum.** For 5 mL. Lyophilized bovine serum containing component concentrations suitable for the quality control of the clinical laboratories, and without preservatives which might interfere with the tests.

### PREPARATION AND USE

1. Open a vial very carefully, avoiding any loss of the lyophilized material.
2. Pipette exactly 5.00 mL of distilled water into the vial. The component values depend on the accuracy of this reconstitution step.
3. Close the vial with the stopper and let stand for 20 minutes at room temperature.
4. Swirl gently, avoiding the formation of foam, to ensure complete dissolution of contents.
5. If the material is to be used for analysis of trace elements, avoid contact of the reconstituted material with the stopper to prevent a possible contamination.
6. The reconstituted control serum is to be treated like the patient serum.

### STORAGE AND STABILITY

Store at 2-8°C.

The lyophilized Biochemistry Control Serum is stable until the expiration date given in the label.

The components of the reconstituted material are stable for at least 7 days at 2-8°C and 30 days at -20°C (when frozen once), excepting:

- Alkaline phosphatase is stable 5 hours at 2-8°C and 30 days at -20°C. It is recommended that the reconstituted material is allowed to stand for 1 hour at room temperature before measurement.
- NEFAs are stable 2 days at 2-8°C and 30 days at -20°C.

CK and bilirubin are sensitive to light. Store the vials protected from light.

Discard the vial if there are signs of microbial contamination or excessive turbidity in the reconstituted product.

### WARNINGS AND PRECAUTIONS

Exercise the normal precautions required for handling all laboratory reagents. Safety data sheet available for professional user on request. Disposal of all waste material should be in accordance with local guidelines. Any serious incident that might occur in relation to the device shall be reported to BioSystems S.A.

### SPECIFIED VALUES

The assigned concentration values for components are batch-dependent and their traceability is shown in the enclosed value sheets. Traceability of the results can be assured only if the BioSystems reagents and recommended measurement procedures are used.

The suggested intervals of acceptable values have been calculated from previous experience in interlaboratory variability and are given for orientation only. Each laboratory should establish its own precision parameters.

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